

**MICROTITER MICROPLATE SURFACE SELECTION GUIDE**

<b>Detection Method</b>	<b>Cat#</b>	<b>Description</b>	<b>Binding</b>	<b>Suggested Molecules to Bind</b>
		<b>Medium Binding Clear Plates and Strips</b>		
<b>Colorimetric</b>	9502227	96 Well Plate UB, RE bottom	Hydrophobic	molecules with large hydrophobic regions, polyclonal antibodies, proteins, glycolipids, ovalbumin lipopolysaccharides, lipoproteins, immunoglobulins, whole cell preparations, Staphylococcus protein A
	3355	Immulon 1B Plate, flat bottom		
	3555	Immulon 1B Plate, round bottom		
	95029350	Strip 1x8 assembled UB, RE-bottom		
	95029100	Strip 1x8 assembled EB, RE-bottom		
	6310	Immulon 1 B 1x12 Strip assembled, flat bottom		
	6301	Immulon 1 B 1x12 Strip, flat bottom		
	6505	Immulon 1 B 2x8 Strip assembled, flat bottom		
95029390	Breakable Strip 1x8 Assembled UB, RE-bottom			
		<b>High Binding Clear Plates and Strips</b>		
<b>Colorimetric</b>	95029330	96 Well Plate EB, RE bottom	Hydrophilic	positively charged molecules, proteins, peptides, glycopeptides, glycoproteins, polyclonal and monoclonal antibodies, immunoglobulins, purified antigens, haptens, whole cell preparations
	3455	Immulon 2 HB Plate, flat bottom		
	3655	Immulon 2 HB Plate, round bottom		
	3855	Immulon 4 HBX Plate, flat bottom		
	6309	Immulon 2 HB 1x12 Strip assembled, flat bottom		
	6302	Immulon 2 HB 1x12 Strip, flat bottom		
	6506	Immulon 2 HB 2x8 Strip assembled, flat bottom		
	6405	Immulon 4 HBX 1x12 Strip assembled, flat bottom		
	6404	Immulon 4 HBX 1x12 Strip, flat bottom		
	6508	Immulon 4 HBX 2x8 Strip assembled, flat bottom		
95029180	Breakable Strip 1x8 Assembled EB, RE-bottom			
		<b>Medium Binding White Plates and Strips</b>		
<b>Luminometric Fluorometric</b>	9502887	White 96 Well Plate UB, RE bottom	Hydrophobic	molecules with large hydrophobic regions, polyclonal antibodies, proteins, glycolipids, ovalbumin lipopolysaccharides, lipoproteins, immunoglobulins, whole cell preparations, Staphylococcus protein A
	7416	Microlite 1 Plate, flat bottom		
	7571	Microlite 1+ Plate, flat bottom		
	7705	Microfluor 1 White Plate, flat bottom		
	6905	Microfluor 1 White Plate, round bottom		
	95029510	White Strip 1x8 assembled UB, RE-bottom		
	95029660	White Breakable Strip 1x8 Assembled UB, RE-bottom		
	7421	Microlite 1 1x12 Strip assembled, flat bottom		
	7403	Microlite 1 1x12 Strip, flat bottom		
	7561	Microlite 1+ 1x12 Strip assembled, flat bottom		
7566	Microlite 1+ 1x12 Strip, flat bottom			
		<b>High Binding White Plates and Strips</b>		
<b>Luminometric Fluorometric</b>	7417	Microlite 2 Plate, flat bottom	Hydrophilic	positively charged molecules, proteins, peptides, glycopeptides, glycoproteins, polyclonal and monoclonal antibodies, immunoglobulins, purified antigens, haptens,
	7572	Microlite 2+ Plate, flat bottom		
	7905	Microfluor 2 White Plate, flat bottom		
	7105	Microfluor 2 White Plate, round bottom		
	7400	Microlite 2 1x12 Strip assembled, flat bottom		
	7410	Microlite 2 1x12 Strip, flat bottom		
	7562	Microlite 2+ 1x12 Strip assembled, flat bottom		
	7567	Microlite 2+ 1x12 Strip, flat bottom		
		<b>Medium Binding Black Plates and Strips</b>		
<b>Fluorometric</b>	9502867	Black 96 Well Plate UB, RE bottom	Hydrophobic	molecules with large hydrophobic regions, polyclonal antibodies, proteins, glycolipids, ovalbumin lipopolysaccharides, lipoproteins, immunoglobulins, whole cell preparations, Staphylococcus protein A
	7605	Microfluor 1 Black Plate, flat bottom		
	7005	Microfluor 1 Black Plate, round bottom		
	95029450	Black Strip 1x8 assembled UB, RE-bottom		
		<b>High Binding Black Plates</b>		
<b>Fluorometric</b>	7805	Microfluor 2 Black Plate, flat bottom	Hydrophilic	positively charged molecules, proteins, peptides, glycopeptides, glycoproteins, polyclonal and monoclonal antibodies, immunoglobulins, purified antigens, haptens, whole cell preparations
	7205	Microfluor 2 Black Plate, round bottom		
		<b>Streptavidin Coated Plates</b>		
<b>Colorimetric</b>	95029263	BioBind Strip 1x8 Assembled	Streptavidin-Biotin bond	biotinylated molecules
	95029293	BioBind Breakable Strip 1x8 Assembled		
<b>Luminometric</b>	95029273	White BioBind Strip 1x8 Assembled		
	95029303	White BioBind Br. Strip 1x8 Assembled		
<b>Fluorometric</b>	95029283	Black BioBind Strip 1x8 Assembled		